

PONDICHERRY UNIVERSITY

(A Central university)



SCHOOL OF ENGINEERING AND TECHNOLOGY

DEPARTMENT OF COMPUTER SCIENCE

Master of Computer Application

NAME : MUHAMMED BINSHAD P T

REG. NO. : 23370038

SEMESTER : THIRD SEMESTER

SUBJECT : ISM

**TITLE : COMPREHENSIVE RISK ASSESSMENT AND
MITIGATION FOR IT ASSETS IN UNIVERSITY CHEMISTRY
LAB**

INDEX

S.no	Asset
1	Introduction
2	Computers And Laptops
3	Electronic Scales (Digital Balances)
4	Spectrophotometers
5	Data Loggers
6	Laboratory Management Software
7	Fume Hood Controllers
8	Printers
9	Data Storage Devices (External Hard Drives, Cloud Storage)
10	Electronic Pipettes
11	UPS Systems (Uninterruptible Power Supplies)
12	Barcode Scanners

INTRODUCTION

In today's academic landscape, IT assets are fundamental to advancing research and enhancing operational efficiency within university labs. Chemistry labs, in particular, rely heavily on technology to conduct experiments, analyze data, and manage sensitive information. However, with the integration of these tools comes the potential for risks that can disrupt operations, compromise data integrity, or pose safety hazards. This document provides a comprehensive assessment of 11 essential IT assets found in Pondicherry university Chemistry lab, detailing their roles, associated risks, and mitigation strategies. Implementing these measures ensures the secure, reliable, and efficient functioning of lab environments.

1. Computers and Laptops

Computers and laptops are essential to the lab's digital infrastructure, enabling researchers to conduct simulations,

perform data analysis, maintain records, and access critical chemical databases. They store and process vast amounts of sensitive data, making them indispensable in managing complex experiments and producing accurate results.

Name: Dell Precision 5560

Specifications: Intel Core i7-11800H, 16GB RAM, 512GB SSD, NVIDIA Quadro T1200, Windows 10 Pro, 15.6" Full HD display, and USB 3.2, Thunderbolt 4, and HDMI ports.

- **Risk Assessment:** These devices face significant risks such as malware attacks, unauthorized access, and accidental hardware failures, which could lead to data loss and impact research continuity.
- **Risk Mitigation:** Installing antivirus software and firewalls can block malicious threats, while regular backups to secure off-site or cloud storage protect against data loss. Multi-factor authentication (MFA) and strong passwords further secure access to sensitive data, safeguarding research integrity.

2. Electronic Scales (Digital Balances)

Electronic scales provide highly accurate measurements of chemicals and substances, which is essential for conducting precise and replicable experiments. They are sensitive to

calibration errors and physical damage, which can compromise measurement accuracy.

Name: Sartorius Entris II Essential Analytical Balance

Specifications: Capacity of 220g, readability of 0.1mg, built-in calibration, and LCD backlit display with touchscreen operation. Interface includes RS-232 and USB for data transfer.

- **Risk Assessment:** Electronic scales are vulnerable to calibration errors, physical damage, and power fluctuations, which can lead to inaccurate measurements and affect experiment outcomes.
- **Risk Mitigation:** Regular calibration and maintenance help maintain accuracy, while protective covers safeguard the equipment from spills. Using stable, surge-protected power sources helps prevent damage from power surges, ensuring the scales perform reliably.

3. Spectrophotometers

Spectrophotometers analyze the light absorption and emission properties of chemical samples, allowing for in-depth analysis critical to chemical research. This sophisticated equipment

requires regular maintenance to ensure reliable performance and accuracy.

Name: Thermo Scientific GENESYS 180 UV-Vis Spectrophotometer

Specifications: Wavelength range of 190-1100 nm, dual-beam optics, touchscreen control, 1.8 nm spectral bandwidth, and USB and Ethernet connectivity.

- **Risk Assessment:** Spectrophotometers are prone to calibration errors, software malfunctions, and physical degradation, all of which can lead to data inaccuracies and compromise experiment validity.
- **Risk Mitigation:** Routine preventive maintenance and calibration checks enhance accuracy, while regular software updates reduce risks from bugs. Ensuring proper storage and handling also prevents damage, and training users on best practices minimizes the likelihood of misuse.

4. Data Loggers

Data loggers record environmental conditions like temperature and pressure, providing essential data for experiments requiring consistent monitoring. These devices are integral to

maintaining the stability of certain experimental parameters over time.

Name: HOBO UX100-011 Temperature Data Logger

Specifications: Temperature range -20°C to 70°C, accuracy of $\pm 0.21^\circ\text{C}$, 84,650 data points, USB data offload, and LCD screen. It also has a battery life of up to 1 year.

- **Risk Assessment:** Data loggers face risks like data corruption, hardware malfunctions, and software glitches, potentially leading to the loss of critical data or inaccurate recordings.
- **Risk Mitigation:** Storing data redundantly in multiple locations, including secure cloud storage, reduces the risk of data loss. Encryption protects sensitive data, and regular exports and backups ensure accessibility. Routine functionality checks further help in identifying issues early.

5. Laboratory Management Software

Laboratory management software facilitates tracking experiments, managing inventory, and storing confidential information. It centralizes lab management tasks and supports an organized workflow, making it critical for lab efficiency.

Name: LabWare LIMS (Laboratory Information Management System)

Specifications: Cloud-based, customizable modules for sample management, test scheduling, data tracking, and report generation. Supports multi-user environments, security features, and integration with lab instruments

- **Risk Assessment:** Laboratory management software is vulnerable to cybersecurity threats like hacking and data breaches, which can lead to data loss or exposure of sensitive information.
- **Risk Mitigation:** Implementing encryption, role-based access control, and regular security updates helps maintain data security. Conducting vulnerability assessments identifies potential risks, while cybersecurity training for lab personnel promotes safe practices, reducing the likelihood of breaches.

6. Fume Hood Controllers

Fume hood controllers regulate airflow within fume hoods, essential for venting hazardous fumes and protecting lab

personnel. These controllers help maintain a safe environment by preventing harmful exposures.

Name: AirClean Systems Fume Hood Controller

Specifications: Digital control for airflow adjustments, LCD touch screen, temperature and humidity sensors, alarms for airflow fluctuations, and connectivity options for data monitoring.

- **Risk Assessment:** Malfunctions in fume hood controllers could expose staff to toxic substances, posing serious health risks if the hood fails to operate correctly.
- **Risk Mitigation:** Regular maintenance and inspections ensure that controllers function correctly. Installing alarms and fail-safes provides immediate notification of system failures, allowing prompt action. Securing access to these controls also prevents unauthorized tampering, enhancing lab safety.

7. Printers

Printers are used to produce hard copies of lab reports, data, and experimental documentation, supporting physical record-

keeping and data sharing. They ensure that vital documentation is available and easily accessible.

Name: HP LaserJet Pro MFP M428fdw

Specifications: Monochrome laser printer, print speed of up to 40 ppm, 600 x 600 dpi resolution, supports duplex printing, and has Ethernet, Wi-Fi, and USB connectivity.

- **Risk Assessment:** Printers face risks like unauthorized access, which could lead to data exposure, and mechanical issues like paper jams or hardware malfunctions, potentially disrupting lab workflows.
- **Risk Mitigation:** Securing print queues with PIN access controls restricts unauthorized access to sensitive documents. Regular maintenance and software updates prevent common printer issues, and data encryption further protects any sensitive information.

8. Data Storage Devices (External Hard Drives, Cloud Storage)

Data storage devices, such as external hard drives and cloud storage, play a critical role in storing large volumes of lab data,

offering reliable backups and easy access for analysis. They are essential for safeguarding the lab's experimental data.

Name: Seagate Expansion Portable External Hard Drive 2TB

Specifications: 2TB storage, USB 3.0 interface, plug-and-play functionality, and compatible with Windows and Mac.

- **Risk Assessment:** These storage devices are vulnerable to physical damage, theft, and cybersecurity threats, leading to data breaches or loss of valuable research information.
- **Risk Mitigation:** Encrypting all data stored on external drives prevents unauthorized access, while securely storing drives in fireproof locations minimizes physical risks. Using reputable cloud providers with advanced security protocols and employing strong access controls protects data stored offsite.

9. Electronic Pipettes

Electronic pipettes enable precise dispensing of small liquid volumes, ensuring accuracy in experiments that require exact measurements. They are widely used to maintain consistency in sample preparation.

Name: Eppendorf Research Plus Electronic Pipette

Specifications: Volume range of 0.5-10 μL , fully autoclavable, digital display for accuracy, battery-powered with rechargeable capability, and multiple modes for different dispensing needs.

- **Risk Assessment:** Miscalibration, wear, or malfunction in electronic pipettes can lead to inaccurate measurements, which could compromise experimental outcomes.
- **Risk Mitigation:** Regular calibration and inspections maintain pipette accuracy, while following manufacturer handling guidelines reduces the risk of damage. Storing pipettes in a clean, controlled environment when not in use extends their lifespan and helps retain measurement accuracy.

10. UPS Systems (Uninterruptible Power Supplies)

UPS systems provide backup power to critical lab equipment, protecting against data loss and equipment damage during power outages. They help ensure the lab operates continuously, even in the event of power interruptions.

Name: APC Smart-UPS 1500VA

Specifications: 1500VA/1000W, LCD display, pure sine wave output, battery backup with runtime up to 10 minutes at full load, surge protection, and USB and Ethernet connectivity for monitoring.

- **Risk Assessment:** UPS failures could cause data loss or equipment damage, disrupting experiments and potentially leading to inaccurate results.
- **Risk Mitigation:** Regular testing and maintenance ensure UPS reliability, while load management prevents overloading. Routine battery replacements ensure the system remains functional during outages, and having a secondary backup source for essential equipment provides additional security.

11. Barcode Scanners

Barcode scanners streamline inventory management, enabling fast and accurate tracking of chemical containers and samples. These devices enhance efficiency in organizing and managing lab resources.

Name: Zebra DS2208 Barcode Scanner

Specifications: 1D and 2D scanning capabilities, USB interface, plug-and-play setup, ergonomic design, scan speed of up to 100 scans per second, and high-resolution scan accuracy.

- **Risk Assessment:** Scanners are at risk of physical damage, software malfunctions, or malware infection, potentially resulting in inventory errors or data loss.
- **Risk Mitigation:** Regular software updates reduce risks from bugs or vulnerabilities, and physical cases protect scanners from damage. Conducting routine device checks ensures accurate performance, while antivirus software guards against malware that could compromise data accuracy.

References

1. www.sciencedirect.com

2. www.labmanager.com
3. www.thermofisher.com
4. www.bhphotovideo.com
5. www.eppendorf.com
6. www.labwrench.com

